

Parts Catalogue

Therapy Products



Parts

If you need further information regarding the parts for your Laerdal products, please contact us. Trained personnel are available to answer your questions and provide additional technical data.

All parts listed for a specific product might not be available in all markets

*Product specifications are subject to change without notice.
Actual colours on the parts might differ from the colours in this catalogue.*

For online purchase please visit

www.laerdal.com



stretcher

Contents

HeartStart MRx	4-7
HeartStart FR – FR2+	8
HeartStart HS1-FRx	9
HeartStart FR3	10
Laerdal Silicone Resuscitators	12-13
Pocket Mask	14
Paediatric Pocket Masks	14
Stifneck Select	14
SpeedBlocks Head Immobilizer	14
Baxstrap Spineboard	15
V-VAC	15
Laerdal Suction Unit	16-17
Laerdal Compact Suction Unit LCSU4	18



M3541A
Carry case for HeartStart
MRx, incl. 3 accessory
pouches/shoulder strap



989806100851
MRx Ambulance Wall Mount



M3538A
MRx battery



40457C Recording Paper Rolls For Cm100
(Pack Of 10 Rolls)
40457D 50mm printing paper (HS 3000QR,
HS 4000, MRx) (80 rolls)
989803138171 MRx 75mm printer paper 10 rolls
989803138181 MRx 75mm printer paper 80 rolls



948003 Power cable



948002 DC Power
Supply for MRx




M3725A Test load for use with
M3508A cable (Plug connector)
M1781A Test load for use with
M3507A cable (Barrel connector)
94 50 02 HeartStart FR to Symbio Tester
adapter cable (not for patient use)

Temperature probes for HeartStart MRx (M3536A/M3535A):

- 21078A** Reusable Skin Surface Temperature Probe
- 21082A** Reusable connector cable for disposable
temperature probes.
- 21090A** Disposable Oesophageal/Rectal Temperature Probe
1 box = 20 probes. (Requires 21082A)
- 21091A** Disposable Skin Surface Temperature Probe
1 box = 20 probes. (Requires 21082A)



ECG Cables

ECG CABLES		ECG 5 LEAD ECG	
	989803147691 Trunk Cable (1,3 meters) (IEC) (standard from Q3/06) <i>The old long trunk cable (2,7m) can be ordered as sparepart (M1949A)</i>		M1971A 5-lead ECG patient cable (extremities. + V1) M1978A 5-lead ECG patient cable – chest (V2 to V6) (To be used together with M1971A)
	M3526A 3-lead ECG patient cable (snap) (1,6 meter) (IEC)		
	M3529A 7-lead ECG patient cable (snap) (IEC) (To be used together with M3526A)		

SpO₂

M1191A SpO ₂ Adult (>50 kg) 	M1194A SpO ₂ Ear clip (>50 kg) 	Nelcor Disposable SpO₂ Sensor  M1943A SpO ₂ adapter to Nelcor single use probe (M1902/4B) (1,1 meter) M1904B Adult (>50kg) M1903B Pediatric (15-50kg) M1902B infant (4-15kg)
 M1192A SpO ₂ Child (15-50 kg)	 M1195A SpO ₂ Infant (4-15 kg)	
 M1941A SpO ₂ Adapter extension cable (2 meter)m		

NBP Antimicrobial Reusable Cuffs

	
M4552B Easy Care Cuff, Infant (10 - 15 cm Circ.) M4553B Easy Care Cuff, Pediatric (14 - 21.5 cm Circ.) M4554B Easy Care Cuff, Sm Adult (20.5 - 28.5 cm Circ.) M4555B Easy Care Cuff, Adult (27.5 - 36 cm Circ.) M4556B Easy Care Cuff, Adult XL (27.5 - 36 cm Circ.) M4557B Easy Care Cuff, Large Adult (35 - 45 cm Circ.) Cuff cover length 66.5 cm. M4558B Easy Care Cuff, Large Adult XL (35 - 45 cm Circ.) Cuff cover length 83 cm	M1598B 1,5 meter NIBP tubing M1599B 3 meter NIBP tubing

Reusable Comfort Cuffs

	M1572A Pediatric size cuff M1573A Adult small size cuff M1574A Adult size cuff M1575A Large Adult size cuff M1576A Thigh size cuff M1578A Multi-Patient Reusable Cuff Kit (M1573A, M1574A, M1575A, M1576A)
--	--


HeartStart MRx Cat.No. M3536A

External Defibrillation Paddles



M3543A
Water resistant external anterior/
anterior paddle set with patient
contact indicator (PCI)

Pads Cable



M3508A
Pads Cable – (Plug connector)



M3507A
Pads Cable – (Barrel connector)



M4763A
Therapy/CPR Cable

Adapters

05-10200
Pads Adapter
(Barrel connector to plug style)

05-10100
Zoll™ Pads Adapter
Can be used with Zoll 1200, 1400, 1600
and M-Series



05-10000
QUIK-COMBO Pads Adapter:
Can be used with Medtronic Physio-
Control LifePak 9, 10C, 11, 12.



989803146981
256MB MRx Data
Card and Tray



M3854A
Datacard and tray



Kat. No	Description	CO ₂ Only	O ₂ /CO ₂	Intubated/ LMA or Combitube	Nasal Only	Oral/Nasal	Duration for use	Primary Applications
M2520A	Pediatric Smart CapnoLine™ O ₂		•			•	8H	10-20 kg
M2522A	Adult Smart CapnoLine™ Plus O ₂		•			•	8H	>55 kg
M2524A	ETCo2 Pediatric Smart Capnoline Single purpose	•				•	8H	10-20 kg
M2526A	ETCo2 Adult Smart Capnoline Single purpose	•				•	8H	>55 kg
M1920A	Adult/Pediatric FilterLine® Set	•		•			6H	>2 kg
M1923A	Infant/Neonatal	•		•			72H	>2 kg

Dual purpose = a probe that enable monitoring of CO₂, and enable supply of O₂.
Humidified = can be used up to 24 hours. All others can be used to 8 hours.



M2520A
Peadiatric Smart Capnoline O₂

M2522A
Adult Smat Capnoline O₂



M2526A
Adult Smart Capnoline




M1920A
Adult/Pediatric MicroStream
Filterline

Therapy (Shock & Pacemaker) electrodes

Electrodes for FR2, MRx, HS 4000



989803158211 Adlt Defibrillation Pads (1 Pair)
989803158221 Adult Defibrillation Pads (5 Pairs)
M3713A
Adult/Child Plus with plug style connector



M3717A
Infant Plus with plug
style connector

Electrodes for MRx, HS 4000
(with barrel style connector)



M3502A Adult Multi Elect For Use With
Cm100/Heartstart4000 (Pack Of 10)



9851-008
Symbio CS1201 (12-lead rhythm
simulator tester)

9801-008
Symbio CS301 (3-lead rhythm
simulator tester)



989803162401
MRx Replace CPR Meter



801-10850 Patient
Adhesives (10x3)

NOT ILLUSTRATED

989806100852 Lower Mount "A" Bracket	989803135341 4 Bay Battery support system for 2 sealed lead acid and 2 lithium ion batteries
989806100853 Lower Mount "B" Bracket	861440 HeartStart Telemedcare System Classic Care Edition (English Only)
989806100854 Upper Mount, Top Bracket	989803146981 256 MB MRx Data card and Tray
989803163291 Patient Adhesive (10pc.) DP2	989803158081 5-lead set, limb, IEC
M1191B Reusable Adult SpO2 finger sensor	989803158091 5-lead set, chest, IEC
M3716A Adult/Child Radiolucent with plug style connector	989803158661 Repl. Pads/CPR Cable
M3718A Adult/Child Radiotransparent (Reduced Skin Irritation) with plug style connector	989803172501 MRx Black Soft Carry Bag
M3719A Infant Radiotransparent (Reduced Skin Irritation) with plug style connector	989803176161 5 lead set, snap, shield AAMI, Chest
21075A Esophageal/rectal temp probe, Adult Reusable	989803176171 5 lead set, snap, shield AAMI, Limb
M3527A 7-lead chest lead set, snap, AAMI 321040	989803176181 5 lead limb lead wires
989803135301 2-BAY ANALYZER FOR HS Li-ION BATTERIES	989803176191 5 lead chest lead wires
989803135331 Battery charger: 4 Bay Battery support system for 4 lithium ion batteries	M1193A Neonatal Hand/Foot SP02 Sensor
	M1530A 5 Lead ECG Cable, IEC
	M1663A 10-ELECTRODE CABLE SETS - 5 WIRE
	M2202A Adult radiotranslucent foam solid gel electrode Size: 2" (50mm) diameter
	M3501A Adult Pads with barrel connector
	M3504A Infant Pads with barrel connector
	M3539A MRx Power Module
	M3542A External Paddles
	M3545A MRx Data Card and Tray
	M5521A MRx handle (Green)
	M5522A MRx handle (Blue)
	M5525A MRx Handle (grey)
	M5529A Heartstart MRx DC Power Module
	M1510A 3 lead ECG Patient Trunk Cable, IEC
	M1615A Shielded 3-Lead Set, Snaps, IEC
	M1633A Shielded 5-Lead Set, Grabbers, IEC
	M1635A Shielded 5-Lead Set WITH Safety Conn.
	M1669A 3 LEAD TRUNK CABLE AAMI/IEC 2.7m
	M3528A 3-electrode limb lead set, snap, IEC
	M3549A Wide Bed Rail Hook (Metal)
	M4737A Keyboard Overlay MRx
	M5530A Combiner plug for 3-wire lead set
	989803166021 Adult / Child Pre-connect
	989803158061 5-lead set, limb, AAMI
	989803158071 5-lead set, chest, AAMI
	M4612A Adult solid get ECG electrode Shape: teardrop
	M4613A Adult solid get ECG electrode Shape: teardrop
	M4761A Q-CPR Compression Sensor
	M4762A CPR Assistent Replacement Pads (10)
	989803185181 Black Soft Carry Case
	989803176541 HeartStart MRx Quick Disconnect DC Power Cable

HeartStart FR - FR2+ Cat. Nos. M3841A, M3840A


Batteries



M3863A
Long-Life LiMnO2
Battery. Disposable
and recyclable. For
FR2 and FR2+ only.



M3854A
Datacard and tray



M3870A Infant/Child
Reduced-Energy
Defibrillator Pads



07-10900
Training pads



40457C
Recording Paper Rolls
(Pack Of 10 Rolls)



989803158211 Adult Pads (1 set)

989803158221 Adult Pads
(1 case = 5 sets)

M3713A
Multifunction Electrode Pads
Adult Plus (1 case = 10 sets)



M3864A
Training & Administration Pack
for FR2 series Defibrillator



NOT ILLUSTRATED

945021	Label status AED Trainer 2 (pkg. 10)
945025	Board AED Trainer 2
945026	Speaker with cable
945027	Cable AED Trainer 2
945055	Programming kit for AED Trainer 2 (optional software and cable)
945060	Technical manual
DP2	Defibrillation electrodes (2 pair)
M3871A	FR2 Paediatric Training Pads
68-PCHAT	Fast Response Kit
941308	Softpack, FR2 Army Green
941320	FR2 Rescue Rucksack
941404	Compact Flash Card Reader USB
M3853A	Data Card Tray only
M3869A	Vinyl Carrying Case for FR2 AED



941311
Carry Case (Hard Pelican)

Adapters

05-10200
Pads Adapter
(Barrel connector to plug style)

05-10100
Zoll™ Pads Adapter
Can be used with Zoll 1200, 1400, 1600 and M-Series

05-10000
QUIK-COMBO Pads Adapter.
Can be used with Medtronic Physio-Control LifePak 9, 10C, 11, 12.

HeartStart HS1 Cat. Nos. M5066A, M5067A

Shock Electrodes



M5071A
Adult SMART Pads
Cartridge (1 Set)



M5072A
Infant/Child SMART
PADS Cartridge (1 set)



M5073A
Adult Training Pads
with Replaceable Seal

Battery

M5070A
(For HS1 and FRx) This battery is a disposable, lithium manganese dioxide, long-life battery that will typically last for four years with the defibrillator in standby mode.



NOT ILLUSTRATED

M5075A	Infant/child training pads cartridge
M5093A	Replacement adult training pads
M5094A	Replacement infant/child training pads
M5089A	External manikin adapters (10 on a roll)
941320	Rescue Rucksack
941343	HeartStart Wall Box with Alarm
941345	Heartstart kit Barrier, wipe & gloves
989803139531	Carry case
989803143051	Case Capture palm SW
989803170911	AED Green Awareness Placard
989803170931	AED Wall Sign, Green
M3859A	Secure Pull Seal Wall Mount
M5075A	HeartStart HS1 Carry Case
M5076A	HeartStart HS1 Slim Carry Case
M5090A	Adult Pad Placement Guide

HeartStart FRx Cat. Nos. M861304, 861305

Carry Case

941360
FRx carrying case made of wipeable urethane.



Battery

M5070A
(For HS1 and FRx) This battery is a disposable, lithium manganese dioxide, long-life battery that will typically last for four years with the defibrillator in standby mode.



NOT ILLUSTRATED

940314	Bracket without Lock
941345	Heartstart Kit Barrier, wipe & gloves
989803139251	HeartStart FRx Carry Case
989803139281	HeartStart Pad Placement Guide, Infant / Child
989803170891	Wall Mount Bracket
989803170911	AED Awareness Placard (Green)
989803170931	AED Wall Sign, Green
M3858A	AED Wall sign
M3859A	Secure Pull Seal Wall Mount 10 pack
M5090A	Adult Pad Placement Guide

Shock Electrodes



989803139261
SMART Pads II




989803139311
Infant/child key used in infant/child defibrillation with the HeartStart FRx

989803139271
Training Pads II (1 pair)

989803139291
Replacement Training Pads TL (1 pair)





989803150191 FR3 Training Pack



945006 Cable AED to Manikin FR/FR2/FR3



989803149981 SMART Pads III, 1 set
989803149991 SMART Pads III, 5 set



801-10850 Patient Adhesives (10 x 3)



989803150111 Fast Response Kit



989803149941 Q-CPR Upgrade

NOT ILLUSTRATED

- 989803150181 Replacement Pads HS Training Pads III
- 989803150201 Interconnect Cable HS Training Pads III
- 989803149951 Cable link,FR3 to CPR meter
- 989803150041 3-Lead ECG cable,FR3 AAMI
- 989803150051 3-Lead ECG Cable,FR3 IEC
- 989803150241 FR3 Rechargeable Battery For Clinical Use
- 989803157531 Laerdal System Case, Rigid
- 989803162231 CPR meter cradle
- 989803179171 Laerdal System Case Soft, FR3 (W/O Auto-on)
- 989803179191 Laerdal small Soft Case FR3 (W/O Auto-on)
- 989803184831 Laerdal small Soft Case FR3 (W/O Auto-on)
- 989803150241 FR3 Rechargeable Battery For Clinical Use

Laerdal Silicone Resuscitators Cat. Nos. 870050/51/52/53/55,

560200 Non-Rebreathing Valve lv



540105
Disk membranes
(pkg. 10)

540103
Lip valve



870150
Adult Bag 1600MI

851250 Patient valve w/press. reg.



540105
Disk membranes
(pkg. 10)

540103
Lip valve

851252
Pressure Relief Valve 35Cm H2O



860150
Paediatric Ventilation bag,
500 ml

851350 Patient valve w/press. reg. and lock clip



540105
Disk membranes
(pkg. 10)

540103
Lip valve

851252
Pressure Relief Valve 35Cm H2O

851103
Lock clips
(pkg. 10)



850150
Preterm Bag 240MI



530400 Airways, set of 4

TUBE AND CONNECTORS



871000 Silicone extension
tube (28 cm)



850900
Manometer
connector



850500
Expiration
diverter



875400
Intake/Reservoir
Valve



551901 O₂ reservoir,
0.6 litre




531901 O₂ reservoir,
2.6 litres

NOT ILLUSTRATED

54010733	High Efficiency Filter for LSR/TheBag(IE)(25pk)
531906	02 Reservoir 2.6Ltr (Pack Of 50)
551906	02 Reservoir 0.6Ltr (Pack Of 50)
85015005	Preterm ventilation bag, 240ml (J)
86015005	Paediatric ventilation bag, 500ml (J)
87015005	Adult ventilation bag, 1600ml (J)
87540000	Intake Reservoir Valve (Nordic)
87540033	Intake Reservoir Valve (IE)
87540040	Intake Reservoir Valve (EUR)



860300 Display Case Paediatric
870600 Display Case Complete Ad.Mod.
850700 Display Case Preterm



875400 Intake/Reservoir Valve



510103
Cap for LSR Intake
Valve (pkg. 3)



510404
Intake
membranes
(pkg. 10)



871950 Flap valves
(pkg. 2)



531907
Intake Valve Adapter
(23 mm OD)



870120
Hanging loop



511700
Wall bracket



871300
Attachment ring f/
standard head strap



870400 Head strap
w/attachm. ring

Pocket Mask Cat. No. Paediatric Pocket Mask Cat. No.

NOT ILLUSTRATED	
820410	One Way Valve Single
820411	Valve Pocket Mask 100/pkg.
820700	Update kit Pocket Mask
830110	Head Strap P/Mask (Pack Of 10)
820610	Pocket Mask Mask Filter Single
820611	Pocket Mask Mask Filter 10 Pack

Stifneck Select Collars Adult & Cat. No. Stifneck Extrication Collars Cat. No.

NOT ILLUSTRATED	
980700	Stifneck Carrying Bag

SpeedBlocks Head Immobilizer Cat. No.



982100 PadPack
Alignment Pads (Pack 10)

NOT ILLUSTRATED	
983060	SpeedBlocks Universal Base (Qty.1)
983092	SpeedBlocks Block Set (Qty.1)
983096	SpeedBlocks Strap and Pad Replacement Set

BaxStrap Spineboard Cat. No.



982100 PadPack
Alignment Pads (Pack 10)



982512 Premium Strap
w/sewn loop ends



982513
Premium Strap w/
speed clip ends



982517 Spider Straps
for BaxStrap

NOT ILLUSTRATED

- | | |
|-----------|---|
| 701-00040 | Laerdal Pedi-Pad
Backboard Pad (Int) 6pk |
| 980800 | Fullcover Carry Bag
Baxstrap Spineboard |

V-VAC Cat. No.



NOT ILLUSTRATED

- | | |
|--------|---------------------------------|
| 886010 | 300ml Disposable |
| 985001 | V-Vac Replace Cartridge |
| 985002 | Adaptor Tips
(Pack Of 4) |
| 985003 | Male Connectors
(Pack Of 10) |
| 985004 | Catheters (Pack Of 4) |
| 985101 | V-Vac Training Cartridge |
| 985200 | V-Vac Carry Bag
(Yellow) |
| 985030 | V-VAC Handle
Spare Part |

Laerdal Suction Unit Cat. Nos. 780000, 780020, 780030

Reusable Canister System Cat. No. 780000



784000 Reusable canister




Serres Canister System Cat. No. 780030



Serres Canister System Pre 2014



BEMIS canister system (US Only) Cat. No. 780020



Parts for all units Cat. Nos. 780000, 780010, 780020



780433 Strap for Tubing



780434 Hook for Tubing Strap



780436 Fasten bracket left/right




780800
LSU NiMH Battery

CANISTER HOLDER



780432 Release arm,
canister holder



780420 Bemis Canister Holder



780451
Serres Canister Holder



780430
LSU Reusable Canister Holder

NOT ILLUSTRATED


780437	LSU Fastening Bracket
780100	Battery Cable
780120	Power Cable
780130	Carrying handle
780140	Main Rotary Switch
780150	Gasket
780160xx	Cabinet cpl.
780170	Vacuum Connector
780180	Foam Absorbers LSU
780190	Battery Lid
780240	Power Board
780250	MMI Board
780270	Battery Contact w/spring
780500	Motor & Pump Assembly
780510	Pump wear and tear parts
780520	Piston-arm exchange tool
780530	Umbrella valves
780600	LSU mmHg Display
780700	LSU KPa display
780110	Technical manual

Reusable Canister System Pre 20xx

781007	Vacuum
780431	Vacuum canister tubing/ filter tubing 30 cm LSU
781004	Aerosol filter (12 pk)
781005	Aerosol filter (50 pk)

Parts for authorized service personnel only
For other service parts see service manual

Serres Conversion Kit



780931xx Serres Conversion Kit

Reusable Conversion Kit



780900xx LSU conversion kit
reusable canister

Serres Alteration Kit



780930xx Serres Alteration Kit

780210
AC power cord US



780220
AC power cord EU

780230
AC power cord UK

780200
DC power cord



782600 Wall Bracket
w/DC-Power Cord

782610 Wall Bracket
w/AC Power Cord US

782620 Wall Bracket
w/AC-Power Cord EU

782630 Wall Bracket
w/AC-Power Cord UK

782640 Wall Bracket
w/o Power Cord

791900
Socket adapter



Laerdal Compact Suction Unit LCSU4



886104
800 ml Disposable Canister
w/o tubing (Qty: 6)



886102
800 ml Disposable Canister
w/tubing



886100
300 ml Disposable Canister
w/tubing

NOT ILLUSTRATED	
881151	LCSU4, Service Replacement Unit
881152	LCSU4, Service Replacement Unit



886121 Top Cover, w/ Label for control panel

886221 Top Cover, RTCA, w/ Label for control panel.



886105
Patient Tube 1.8 m (6')
Disposable



886106
Vacuum Tube (for 800 ml)



886111
External Battery Charger



886113
Battery, 12V DC NiMH,
Rechargeable



886122 Bottom Cover



884500
12V DC power cord



886128 AC Power-cord UK



886127 AC Power-cord EU



886126 AC Power-cord US



886123
Battery Cover (door)

Compatible Items

Items that can be used with LCSU 2, 3 and 4



			Can be used with:			Version	
Cat. No	Part Name	2nd Part Name	LCSU 4	LCSU 3	LCSU 2	300ml	800ml
886100	300 ml Disp. Canister	w/tubing (Qty.1)	•	•		•	
886102	800 ml Disp. Canister	w/tubing (Qty.1)	•	•	•		•
886104	800 ml Disp. Canister	w/o tubing (Qty.6)	•	•	•		•
886105	Patient Tube 1.8 m (6')	Disposable (Qty.1)	•	•	•	•	•
886106*	Vacuum Tube	(Qty.1)	•	•	•		•
886108	Wire Stand (for 800 ml)		•				•
886111	AC/DC Adapter Charger	w/o Power-cord	•	•	•	•	•
886112	External Battery Charger	w/o Power-cord	•	•		•	•
886113	Battery, 12V DC NiMH	Rechargeable	•	•		•	•
886114	DC Power-cord		•	•	•	•	•
884500	12V DC Power Cord	(1 each)	•	•	•	•	•
881151	LCSU 4, Replacement unit		•			•	•
881152	LCSU 4, RTCA Repl. unit		•			•	•
886121	Top Cover, w/Label for	control panel	•			•	•
886122	Bottom Cover		•			•	•
886123	Battery Cover (door)		•			•	•
886126	AC Power-cord US		•	•	•	•	•
886127	AC Power-cord EU		•	•	•	•	•
886128	AC Power-cord UK		•	•	•	•	•

* Items are not compatible with 1200 ml canister for LCSU 2.



Laerdal

helping save lives